



Date	Est. #				Shift				# of passed tests:		Page		of							
									# of tests taken:											
HEAD DEFECTS	Test #1		Test #2		Test #3		Test #4		Test #5		Test #6		Test #7		Test #8		Test #9		Test #10	
	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.
Improper dressing (hide, horns, eyelids, etc.)																				
Any visible GIT contents on head																				
Parts missing or >50% tongue/l. node/muscle missing																				
Improper preparation of head (e.g. <75% of cheek muscle incised and <3/1 ratio exposed muscle surface to connective tissue, etc.) and palpation of tongues not properly done as per 6.3.2 .																				
Loss of ID																				
Mode/score	TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR	
Accept/reject	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re
Line speed																				
HEAD DEFECTS	Test #11		Test #12		Test #13		Test #14		Test #15		Test #16		Test #17		Test #18		Test #19		Test #20	
	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.	Time	Init.
Improper dressing (hide, horns, eyelids, etc.)																				
Any visible GIT contents on head																				
Parts missing or >50% tongue/l. node/muscle missing																				
Improper preparation of head (<75% of cheek muscle incised and <3/1 ratio exposed muscle surface to connective tissue, etc.) and palpation of tongues not properly done as per 6.3.2 .																				
Loss of ID																				
Mode/score	TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR		TNR	
Accept/reject	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re	Ac	Re
Line speed																				

T = Tightened, N = Normal, R = Reduced, Ac = Accept, Re = Reject